## Q.P. Code : 613300

	(3 Hours)	[ Total Marks	: 80
N.B	<ul> <li>(1) Question No.1 is compulsory.</li> <li>(2) Answer any three questions out of remaining five questions</li> <li>(3) Figure to the right indicate full marks</li> </ul>	ons.	5
(t (c	a) Explain the physiological effects of low frequency currents b) Explain Circus Motion Theory c) How is Bipolar Mode different from Unipolar mode in ESU d) Explain the Principle of Condensor Method		5 5 5 5
	what do you understand by Cardioversion. Explain in detail of a block diagram  Explain with the help of a neat circuit diagram the INST and of Defibrillator Operation	_	10 10
	Explain in detail with the help of a circuit diagram, the working Diathermy Unit, also mention the technical specifications.  Explain in detail the working of a Nerve Muscle Stimulator	g of Shortwave	10 10
	<ul> <li>Explain with the help of a neat block diagram the working of a unit</li> <li>How is temperature monitored and controlled in a Hemodia</li> </ul>		10 10
	Explain with the help of a neat circuit diagram, the CUT BICOAG selections in an ESU machine  A capacitively coupled return electrode dielectric on an ESU1 d=0.0025cm. The dielectric constant of the insulator ε <sub>r</sub> is space ε <sub>0</sub> is 8.8pF/m. Considering the Area as 120cm <sup>2</sup> impedance versus area for the electrode against a person ESU operating at 500kHz	has a thickness 3.2 and free ,compute the	10 10
	) Draw and explain Rate generator circuit and Rate Generator ) Explain with the help of a neat diagram the working o Pacemaker		10 10